

Three-dimensional ultrasound imaging necessitates the use of probes in which, the ultrasound beam can be focused in both the elevation and azimuth direction. This can be attained by using mechanically moving probes matrix arrays or row-column arrays.

Mechanically moving probes use a mechanical translation or rotation of the probe element to cover the volume. This often gives a good image quality for stationary objects, but the probes are mechanically complicated to manufacture and in general are not suitable for fast B-mode and flow imaging.

Traditional matrix array probes [1-5] have $N \times N$ elements, which even for small arrays gives a very large number of transducer connections (1024 for $N=32$). This makes them very difficult to use, as the number of scanner channels required is high (up to four Verasonics Vantage™ 256 systems must be used for one probe), and the data rates are correspondingly large. For example, it also limits their focusing ability, as the full width half max (*FWHM*) resolution is given by $FWHM = l / F\#$, where l is the wavelength, and the F-number ($F\#$) is the ratio between the imaging depth and aperture width. Beam steering necessitates an element pitch of $l/2$ for the matrix array, and perfect focusing for a 32×32 array is, thus, only possible to $16 l$ for a fully populated array for $F\#=1$. Having arrays with the same focusing ability as traditional phased or linear array would increase the element count to $64 \times 64 = 4096$ or even $192 \times 192 = 36,864$, which is impractical in terms of cables and scanner channels.

Row-column arrays (RCA) are a new matrix array type enabling the acquisition of volumetric data with a reduced number of elements [7-17]. An RCA 2-D transducer array consists of two orthogonal 1-D arrays as shown below with connections to either the rows or the columns [7,10,16]. Emissions can be made with either the rows or the columns, and the received channel signal is the combined signals from all elements in either a row or column. The individual matrix elements are not directly addressable, and the number of connections to the probe is therefore $2N$.




Fig. 1: Illustration of a row-column array, which can be interpreted as two orthogonal 1-D arrays. To the right is shown a 2-D transducer array, where each transducer element is addressed by its row or column index, effectively creating the two arrays shown to the left (from [16]).

RCA's can thus be made quite large without a prohibitive number of connections. A 128 x 128 RCA has 256 connections, and the *FHWM* resolution can be maintained down to 64λ , when the element pitch is $\lambda/2$. A 1024 connection RCA would be able to maintain resolution down to 256λ and would have a size 64 times large than a fully populated matrix array with the same number of connections. Such RCA's have a very large penetration depth due to its size, as demonstrated in [18] and near ideal isotropic resolution can be attained with low side lobes [19]. The Point Spread Function (*PSF*) can be symmetric in both the azimuth and elevation direction, with a lateral Full Width Half Maximum (*FHWM*) of approximately 0.6λ , close to the ideal *FHWM* [19].

The imaging strategies developed for other arrays cannot be directly employed on RCA's, as their geometry is different. The elongated elements give edge wave and focusing must be modified. Special beamformers taking this into account have been developed [16, 20] for proper focusing. New imaging schemes have to be devised, as focused emissions for scanning a 3-D volume may be too slow to yield acceptable volume rates. Synthetic aperture schemes have shown to be successful for yielding good B-mode image quality [18] with both a good resolution and contrast [19]. It has also been demonstrated that full velocity imaging can be attained at a high-volume rate with estimating all three velocity components for tensor velocity imaging [21-23], and it is possible to detect low velocity flow with such probes [24-25]. RCA's have also yielded functional imaging in a rat brain using plane wave emission [26]. Volumetric super resolution imaging has also been demonstrated for a 62 + 62 RCA operating at 3 MHz, where the volume resolution was increased by a factor of 18,880 [27].

Row-column arrays are in their early stages of development, where many applications are being investigated for high quality B-mode, velocity, super resolution, and functional imaging. While it is not fully clear what the best imaging schemes are or how the ideal array is fabricated, the arrays show a very significant potential for radically changing volumetric ultrasound imaging, making them very interesting for future research.

Row-column Array Imaging note contributed by Jørgen Jorgen Arendt Jensen from the Center for Fast Ultrasound Imaging, Department of Health Technology at the Technical University of Denmark.

Bibliography

- [1] O. T. von Ramm, S. W. Smith, and H. G. Pavy. High speed ultrasound volumetric imaging system — Part II: Parallel processing and image display. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 38(2):109–115, 1991.
- [2] R. E. Davidsen, J. A. Jensen, and S. W. Smith. Two-dimensional random arrays for real time volumetric imaging. *Ultrason. Imaging*, 16(3):143–163, July 1994.
- [3] S. S. Brunke and G. R. Lockwood. Broad-bandwidth radiation patterns of sparse two-dimensional vernier arrays. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 44(5):1101–1109, Sep 1997.
- [4] J. T. Yen and S. W. Smith. Real-time rectilinear volumetric imaging using a periodic array. *Ultrasound in Medicine & Biology*, 28(7):923–931, Jul 2002.
- [5] A. Austeng and S. Holm. Sparse 2-D arrays for 3-D phased array imaging – Design methods. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 49(8):1073–1086, 2002.
- [7] C. E. Morton and G. R. Lockwood. Theoretical assessment of a crossed electrode 2-D array for 3-D imaging. In *Proc. IEEE Ultrason. Symp.*, pages 968–971, 2003.
- [8] N. M. Daher and J. T. Yen. 2-D array for 3-D ultrasound imaging using synthetic aperture techniques. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 53(5):912–924, 2006.
- [9] S. I. Awad and J. T. Yen. 3-D spatial compounding using a row-column array. *Ultrason. Imaging*, 31(2):120–130, 2009.
- [10] C. H. Seo and J. T. Yen. A 256 x 256 2-D array transducer with row-column addressing for 3-D rectilinear imaging. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 56(4):837–847, April 2009.
- [11] C. E. M. Démoré, A. W. Joyce, K. Wall, and G. R. Lockwood. Real-time volume imaging using a crossed electrode array. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 56(6):1252–1261, 2009.
- [12] A. S. Logan, L. L. P. Wong, and J. T. W. Yeow. 2-D CMUT wafer bonded imaging arrays with a row-column addressing scheme. In *Proc. IEEE Ultrason. Symp.*, pages 984–987, sep 2009.
- [13] A. S. Logan, L. L. P. Wong, A. I. H. Chen, and J. T. W. Yeow. A 32 x 32 element row-column addressed capacitive micromachined ultrasonic transducer. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 58(6):1266–1271, June 2011.
- [14] A. I. H. Chen, L. L. P. Wong, A. S. Logan, and J. T. W. Yeow. A CMUT-based real-time volumetric ultrasound imaging system with row-column addressing. In *Proc. IEEE Ultrason. Symp.*, pages 1755–1758, oct 2011.
- [15] A. Sampaeanu, P. Zhang, A. Kshirsagar, W. Moussa, and R. Zemp. Top-orthogonal-to-bottom-electrode (TOBE) CMUT arrays for 3-D ultrasound imaging. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 61(2):266–276, 2014.
- [16] M. F. Rasmussen, T. L. Christiansen, E. V. Thomsen, and J. A. Jensen. 3-D imaging using row-column-addressed arrays with integrated apodization — Part I: Apodization design and line element beamforming. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 62(5):947–958, 2015.
- [17] T. L. Christiansen, M. F. Rasmussen, J. P. Bagge, L. N. Moesner, J. A. Jensen, and E. V. Thomsen. 3-D imaging using row-column-addressed arrays with integrated apodization — part II: Transducer fabrication and experimental results. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 62(5):959–971, 2015.
- [18] H. Bouzari, M. Engholm, S. I. Nikolov, M. B. Stuart, E. V. Thomsen, and J. A. Jensen. Imaging performance for two row-column arrays. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 66(7):1209–1221, 2019.
- [19] M. Schou, A. S. Havreland, M. Engholm, M. B. Stuart, E. V. Thomsen, and J. A. Jensen. Design of a novel zig-zag 192+192 row column addressed transducer: A simulation study. In *Proc. IEEE Ultrason. Symp.*, pages 1–4, 2018.
- [20] M. B. Stuart, P. M. Jensen, J. T. R. Olsen, A. B. Kristensen, M. Schou, B. Dammann, H. H. B. Sørensen, and J. A. Jensen. Real-time volumetric synthetic aperture software beamforming of row-column probe data. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 68(8):2608–2618, 2021.
- [21] S. Holbek, T. L. Christiansen, M. B. Stuart, C. Beers, E. V. Thomsen, and J. A. Jensen. 3-D vector flow estimation with row-column addressed arrays. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 63(11):1799–1814, 2016.
- [22] M. Schou, L. T. Jørgensen, C. Beers, M. S. Traberg, B. G. Tomov, M. B. Stuart, and J. A. Jensen. Fast 3-D velocity estimation in 4-D using a 62 + 62 row-column addressed array. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 68(3):608–623, 2020.
- [23] L. T. Jørgensen, M. Schou, M. B. Stuart, and J. A. Jensen. Tensor velocity imaging with motion correction. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 68(5):1676–1686, 2020.
- [24] M. Flesch, M. Pernot, J. Provost, G. Ferin, A. Nguyen-Dinh, M. Tanter, and T. Defieux. “4D in vivo ultrafast ultrasound imaging using a row-column addressed matrix and coherently-compounded orthogonal plane waves,” *Phys. Med. Biol.*, vol. 62, no. 11, pp. 4571–4588, 2017.
- [25] J. Sauvage, M. Flesch, G. Ferin, A. Nguyen-Dinh, J. Poree, M. Tanter, M. Pernot, and T. Defieux. “A large aperture row column addressed probe for in vivo 4D ultrafast Doppler ultrasound imaging,” *Phys. Med. Biol.*, vol. 63, pp. 1–12, 2018.
- [26] J. Sauvage, J. Poree, C. Rabut, G. Ferin, M. Flesch, B. Rosinski, A. Nguyen-Dinh, M. Tanter, M. Pernot, and T. Defieux. “4D functional imaging of the rat brain using a large aperture row-column array,” *IEEE Trans. Med. Imag.*, vol. 39, pp. 1884–1893, June 2020.
- [27] J. A. Jensen, M. L. Ommen, S. H. Øygard, M. Schou, T. Sams, M. B. Stuart, C. Beers, E. V. Thomsen, N. B. Larsen, and B. G. Tomov. Three-dimensional super resolution imaging using a row-column array. *IEEE Trans. Ultrason., Ferroelec., Freq. Contr.*, 67(3):538–546, 2020.

Verasonics Inc.

11335 NE 122nd Way, Suite 100, Kirkland, WA 98034

www.verasonics.com | sales@verasonics.com | 425.998.9836



Verasonics reserves the right to change specifications without notice.

© Verasonics, Inc. 2022. All rights reserved. 4002-2022-34