


Next generation multi-channel and phased array ultrasound instrumentation
Engineered for cutting-edge performance in research and industry
32 – 2048 Channels | 20 kHz* – 60 MHz Bandwidth


Advanced Transmit Capability

The transceiver is the heart of an ultrasound acquisition device, but not all transceivers are created equal. Verasonics' proprietary, patented *NXT* transceiver was specifically engineered to provide the next level of transmit and receive fidelity required for cutting-edge ultrasound research and inspection.


- **Industry-leading transmit fidelity:**
Resolution, dynamic range, pulse inversion symmetry, and harmonic distortion.
- **High power transmit:**
Continuous transmit on all channels for milliseconds duration (or seconds depending on configuration).
- **Coded waveform capabilities:**
Advanced waveform coding performance enabled by *NXT* transmit pulse symmetry and repeatability.
- **Benefits of transmit coding:**
 - Increased resolution
 - Higher SNR
 - Faster inspection



No Averages:
64 Acquisitions



64 Averages:
4096 Acquisitions



32-bit Golay Compression:
128 Acquisitions

Waveform Coding Example: 5 MHz TFM Imaging in HDPE, dB Scale

Arbitrary Waveform Generation

- Proprietary algorithms enable generation of high-fidelity arbitrary analog acoustic waveforms using the *NXT* three-level transmitter.
- DAC + linear amplifier level performance with high channel counts for a fraction of the cost and power consumption.



Experimentally measured arbitrary ultrasonic waveform: Pseudo-random 5-11 MHz

Engineered for High-Speed Inspection

Vantage *NXT* features hardware and software engineered for the highest performance in acquisition speed, enabling continuous, uninterrupted data acquisition with scan speeds limited only by acoustic travel time.

- **Ultra-Fast Data Transfer:**
Up to 52.8 Gbit/s (6.6 GB/s) sustained.
- **GPU Direct:**
Data transfer directly to VRAM.
- **Asynchronous Architecture:**
Parallelization of data acquisition, transfer and processing.
- **Large Internal Memory:**
256 MB per Channel.

Vantage *NXT* Specifications

| TRANSMIT SPECIFICATIONS | |
|-------------------------------|---|
| Low Frequency | 20 kHz - 2 MHz* |
| Mid Frequency | 250 kHz - 24 MHz |
| High Frequency | 1.0 MHz - 60 MHz |
| HIFU Operating Band | 250 kHz to 8 MHz |
| Programmable Pulser Amplitude | 3.0 to 192 V p-p |
| Waveform Duration | Standard: ~microseconds Extended Transmit: ~milliseconds HIFU: ~seconds |
| Minimum "on" state duration | 4 ns |
| Resolution (clock frequency) | 2 ns (500 MHz) |
| Voltage transition time | 2-3 ns |
| FET technology | GaN |

| I/O OPTIONS | |
|----------------------------|---|
| Triggers | 4 input/outputs programmable pulse width |
| Quadrature Encoder Inputs* | 3 quadrature encoders |
| Analog Inputs | 2 inputs Signal +/- 5V range Sample rate 4 KHz 12-bit ADC |
| PCIe | Gen 3x8 FireFly™ firefly optical connection with 2 m cable (up to 100 m) |

| RECEIVE SPECIFICATIONS | |
|------------------------------|-------------------------------|
| Low Frequency | 20 kHz - 2 MHz* |
| Mid Frequency | 250 kHz - 60 MHz |
| High Frequency | 250 kHz - 60 MHz |
| ADC sample rate | 125 MHz |
| ADC Resolution | 16-bit |
| Channel memory | 256 MB |
| Sustained data transfer rate | Up to 52.8 Gbits/s (6.6 GB/s) |
| Subsampling | Non-integer, rational factors |
| Data type | 16-bit integer, 16-bit float |

* Future availability

