







L11-5v
Linear Array

Element number: 128 Pitch (mm): 0.3
Elevation focus (mm): 18 Sensitivity (dB): -52 ± 3





L12-3v
Linear Array

Element number: 192 Pitch (mm): 0.2
Elevation focus (mm): 20 Sensitivity (dB): -58.4





C5-2v
Curved Linear Array

Element number: 128 Pitch (mm): 0.506-0.510
Elevation focus (mm): 50-70 Sensitivity (dB): -55.5 -50.5



P4-2v
Phased Array

Element number: 64 Pitch (mm): 0.298-0.302
Elevation focus (mm): 50-70 Sensitivity (dB): -69 -65



NOTE: These transducers have not been cleared by the FDA for diagnostic use.

