




Transducer Specifications L11-5v, C5-2v, P4-2v

Compatible with all Standard Frequency and High Frequency configurations of Vantage systems using the UTA 260-D, UTA 260-S, or the UTA 260-MUX adapters.



Element number: 128 **Pitch (mm):** 0.506–0.510 **Elevation focus (mm):** 50–70 **Sensitivity (dB):** -55.5 – -50.5









C5-2v

Curved Linear Array

Element number: 64 **Pitch (mm):** 0.298–0.302 **Elevation focus (mm):** 50–70 **Sensitivity (dB):** -69 – -65





All technical specifications are provided courtesy of the transducer manufacturer. NOTE: These transducers have not been cleared by the FDA for diagnostic use.



Verasonics reserves the right to change specifications without notice. $\hfill {\ensuremath{\mathbb C}}$ Verasonics, Inc. 2021. All rights reserved. 4000-2020-02